



Excel-based Quantitation Report System (QRS) Trending



■ ***Extension of MassHunter integrated Excel-based reporting***

Excel-based Quantitation Report System (QRS) Videos



Excel-based QRS Trending

Quantitation Report System

- o [Introduction \(Background\)](#) #...
- o [Scenarios \(Background\)](#) #...

Quick Demo

- o [Calibration \(Advanced\)](#) #...
- o [Target Compounds \(Advanced\)](#) #...
- o [Qualifiers \(Advanced\)](#) #...
- o

Advanced Demo

- o [Calibration \(Expert\)](#) #...
- o [Target Compounds \(Expert\)](#) #...
- o [Qualifiers \(Expert\)](#) #...

Technical Design

- o [Design explanation \(Wizard\)](#) #...
- o [Error Recovery Advice \(Wizard\)](#) #...

Videos on
shipping
Quant DVD




Agilent Technologies


MassHunter Quantitative Analysis software
All the performance, all the time

QRS Excel-based Trend Report Value Proposition

- Address trend reporting needs across batches or date/time
 - Monitor samples types
 - Trend by Samples, CALs, QCs, Blanks, Spikes
 - Monitor target compounds or qualifiers results
 - Trend metrics grouped by sample/compound/group
 - Trend organized by time series { before/after, treated/untreated }
- Monitor outliers results
 - Trend by quality metrics



Quantitation Pipeline software architecture




Agilent Technologies

MassHunter Quantitative Analysis
Quantitation Reporting System

Quantitation Pipeline software architecture

Acquisition

Deliverables:
Acquisition data




Quant

Deliverables:
Reporting results




Reporting

Deliverables:
Batch based
Excel-based Reports



RF for QCs over time series




Deliverables:
Trend Reporting
Cross-batch
QA trend metrics




Agilent Technologies

MassHunter Quantitative Analysis
Quantitation Reporting System


Data Warehouse and Analysis Cube (star configuration)



QRS Star Configuration (MultiDimensional Analysis)



Excel-based Quantitation Report System (QRS) Trending



- ***Extension of MassHunter integrated Excel-based reporting***

Welcome to MassHunter Quantitation Report System (QRS).

QRS is an extension of MassHunter Excel-based reporting designed to show trends.
QRS offers specialized ribbons that focus on a particular aspect of the quantitation results.
QRS specialized ribbons offer canned charts that display a wide variety of trends.

QRS offers an interactive stepwise process to building your charts:

1. QRS creates an XML Table from importing multiple XML quantitation reporting results files.
2. QRS creates a PivotTable Cache which organizes the combined cross-batch quantitated results.
3. QRS creates a PivotTable with appropriate measures and filters from the shared PivotTable Cache.
4. QRS creates a PivotTable Chart which displays the trend report graphically.


Select Quantitation Report System Trend Report: _____

Calibration

TargetCompound


TargetQualifier






LevelResponse for ISTD


Looking at particular Calibration levels, are my are my ISTD responses *trending*?






FinalConcentration for SampleType Blank

Looking at particular Target Compounds,
are my Blank concentrations *trending*?







Security Warning Data connections have been disabled Options...

QualifierResponseRatio for SampleType Sample

Looking at particular Qualifier transitions,
are my qualifier ratios *trending*?





Welcome to MassHunter Quantitation Report System (QRS).

QRS is an extension of MassHunter Excel-based reporting designed to show trends.
QRS offers specialized ribbons that focus on a particular aspect of the quantitation results.
QRS specialized ribbons offer canned charts that display a wide variety of trends.

QRS offers an interactive stepwise process to building your charts:

1. QRS creates an XML Table from importing multiple XML quantitation reporting results files.
2. QRS creates a PivotTable Cache which organizes the combined cross-batch quantitated results.
3. QRS creates a PivotTable with appropriate measures and filters from the shared PivotTable Cache.
4. QRS creates a PivotTable Chart which displays the trend report graphically.

Select Quantitation Report System Trend Report: _____

Calibration

TargetCompound

TargetQualifier



Batch	
ProjectID	
BatchID	
SampleID	
AcqDateTime	
AcqMethodFileName	
AcqMethodPathName	
BalanceOverride	
Barcode	
Comment	
Completed	
DADateTime	
DAMethodFileName	

TargetCompound	
ProjectID	
BatchID	
SampleID	
CompoundID	
AccuracyLimitMultiplierLOQ	
AccuracyMaximumPercentDevi...	
AreaCorrectionFactor	
AreaCorrectionMZ	
AverageRelativeRetentionTime	

TargetQualifier	
ProjectID	
BatchID	
SampleID	
CompoundID	
QualifierID	
AreaSum	
CollisionEnergy	
FragmentorVoltage	
GraphicPeakQualifierChromatog...	
IntegrationParameters	
MZ	
MZExtractionWindowUnits	

Calibration trends

Calibration	
ProjectID	
BatchID	
SampleID	
CompoundID	
LevelID	
CalibrationSTDPathN...	
CalibrationType	
LevelAverageCounter	
LevelConcentration	
LevelEnable	
LevelName	
LevelResponse	

Peak	
ProjectID	
BatchID	
SampleID	
CompoundID	
PeakID	
Accuracy	
Area	
AreaCorrectionResponse	
BaselineDraw	
BaselineEnd	
BaselineStandardDeviation	

PeakQualifier	
ProjectID	
BatchID	
SampleID	
CompoundID	
PeakID	
QualifierID	
Area	
BaselineEnd	
BaselineStandardDeviation	
BaselineStart	
Height	



Batch	
ProjectID	
BatchID	
SampleID	
AcqDateTime	
AcqMethodFileName	
AcqMethodPathName	
BalanceOverride	
Barcode	
Comment	
Completed	
DADateTime	
DAMethodFileName	

TargetCompound	
ProjectID	
BatchID	
SampleID	
CompoundID	
AccuracyLimitMultiplierLOQ	
AccuracyMaximumPercentDevi...	
AreaCorrectionFactor	
AreaCorrectionMZ	
AverageRelativeRetentionTime	
AverageResponseFactor	
BlankResponseOffset	

TargetQualifier	
ProjectID	
BatchID	
SampleID	
CompoundID	
Intensity	
MZ	
MZExtractionWindowUnits	

Calibration	
ProjectID	
BatchID	
SampleID	
CompoundID	
LevelID	
CalibrationSTDPathN...	
CalibrationType	
LevelAverageCounter	
LevelConcentration	
LevelEnable	
LevelName	
LevelResponse	

Peak	
ProjectID	
BatchID	
SampleID	
CompoundID	
PeakID	
Accuracy	
Area	
AreaCorrectionResponse	
BaselineDraw	
BaselineEnd	
BaselineStandardDeviation	

PeakQualifier	
ProjectID	
BatchID	
SampleID	
CompoundID	
PeakID	
QualifierID	
Area	
BaselineEnd	
BaselineStandardDeviation	
BaselineStart	
Height	

Target
Compound
trends



Batch	
ProjectID	
BatchID	
SampleID	
AcqDateTime	
AcqMethodFileName	
AcqMethodPathName	
BalanceOverride	
Barcode	
Comment	
Completed	
DADateTime	
DAMethodFileName	

TargetCompound	
ProjectID	
BatchID	
SampleID	
CompoundID	
AreaSum	
AverageRelativeRetentionTime	
AverageResponseFactor	
BlankResponseOffset	

TargetQualifier	
ProjectID	
BatchID	
SampleID	
CompoundID	
QualifierID	
AreaSum	
CollisionEnergy	
FragmentorVoltage	
GraphicPeakQualifierChromatog...	
IntegrationParameters	
MZ	
MZExtractionWindowUnits	

Calibration	
ProjectID	
BatchID	
SampleID	
CompoundID	
LevelID	
CalibrationSTDPathN...	
CalibrationType	
LevelAverageCounter	
LevelConcentration	
LevelEnable	
LevelName	
LevelResponse	

Peak	
ProjectID	
BatchID	
SampleID	
CompoundID	
PeakID	
Accuracy	
Area	
AreaCorrectionResponse	
BaselineDraw	
BaselineEnd	
BaselineStandardDeviation	

PeakQualifier	
ProjectID	
BatchID	
SampleID	
CompoundID	
PeakID	
QualifierID	
Area	
BaselineEnd	
BaselineStandardDeviation	
BaselineStart	
Height	

Qualifier trends





Pivot Table Data
Crunching for
Microsoft Office
Excel 2007
consolidates all the
best functionality of
pivot tables
including Using VBA
to Create Pivot
Tables



Agilent Technologies

QRS_TrendCharts_B_04_00.xlsm - Micro... PivotTable Tools

Home Insert Page Layout Formulas Data QRS Review View Developer Add-Ins Options Design

Identify Batches Create Table Clear All Final Conc. Batch Sample Category

Sample Blank DbIBlk CAL QC CC MatrixSpike MatrixSpikeDup MatrixSpikeBlank

Sample Group Target ISTD Matrix Manually Integrated Outliers Collapse Batch

Create Project Value Category Sample Type Trend Charts Compound Type Charts Outliers

PivotTable Field List

Choose fields to add to report:

- ☒ BatchName
- ☒ SampleName
- ☐ InstrumentName
- ☐ AcqDateTime
- ☐ UserDefined
- ☐ SampleGroup
- ☒ SampleType
- ☐ SamplingDateTime
- ☐ SamplingTime
- ☒ CompoundName
- ☒ Compound...
- ☐ CompoundGroup
- ☐ Integrator
- ☐ IntegrationMetri...
- ☐ Area
- ☐ Height
- ☐ CalculatedConce...
- ☒ FinalConcentra...
- ☐ ManuallyIntegra...

Report Filter

Column Labels

- SampleType
- CompoundName
- CompoundType

Row Labels

- BatchName
- SampleName

Σ Values

- Average of FinalC...

De... Update

Pivot Chart (CRT-)

Pivot Table (PT-)

Pivot Table (PivotTable)

	A	B	C	D
1	Average of FinalConcentration			
2				
3				
4				
5	QRS			
6				
11				
12	QRS-Day4	Sample-1	12.609	11.596
13		Sample-2	0.000	0.000
14		Sample-3	0.000	0.000
15	QRS-Day5	Sample-1	0.000	0.000
16		Sample-2	4.385	4.715
17		Sample-3	0.000	0.000
18		Sample-1	5.127	5.294
19		Sample-2	6.970	6.552

TargetCompound CRT-Sample PT-Sample PivotTable

Quantitation Report System (QRS) : Sample ...

100%

Sharing the Pivot Cache

Many times, you may have to analyze the same dataset in multiple ways. In most cases, this process requires you to create separate pivot tables from the same data source. You might remember that each time you create a pivot table, you are storing a snapshot of your entire dataset in a pivot **cache**. Every pivot **cache** that is created increases your memory usage and file size. Because of this increase in file size, you should consider sharing your pivot **cache**. That is, in those situations in which you need to create multiple pivot tables from the same data source, use the same pivot **cache** to create multiple pivot tables. By using the same pivot **cache** for multiple pivot tables, you gain a level of efficiency when it comes to memory usage and files size.

In previous versions of Excel, when you created a pivot table, Excel would automatically use the same pivot **cache** if it was already being used in another pivot table. In Excel 2007, Excel gives you no choice. In Excel 2007, a new pivot **cache** is created for each pivot table. This means that multiple pivot **caches** exist for the dataset being used. The size of the file bloats with redundant data each time you create a new pivot table.

You can easily work around this potential problem. By simply copying a pivot table and pasting it into a new pivot table, *without* duplicating the pivot **cache**, as you want to the same pivot **cache**, with a different **Qualifier**.

Pivot Cache shared by
Pivot Tables. Different
Pivot Caches created for
Calibration
TargetCompound
Qualifier

Side Effects of Sharing a Pivot Cache

It's important to note that there are a few side effects to sharing a pivot **cache**. For example, suppose you have two pivot tables using the same pivot **cache**. Certain actions affect both pivot tables. They include



Microsoft Visual Basic - QRS_TrendCharts_B_04_00.xlsm - [QRSInitialize (Code)]

File Edit View Insert Format Debug Run Tools Add-Ins Window Help Type a question for help

Ln 22, Col 67

Project - VBAProject

- ThisWorkbook
- Forms
 - QRSTrendReportForm
- Modules
 - GeneratePivotChart
 - QRSCollapseOutline
 - QRSCreatePivot
 - QRSDelete
 - QRSImportXML
 - QRSInitialize**
 - QRSSetMultiDimensions
 - UIAnalysisDelivery
 - UICALCalibrationCurve
 - UICALCalibrationOutlier
 - UICALCalibrationPoints
 - UICALRespFactorTrend
 - UIQualifierOutlier
 - UISHARE
 - UISHARECleanUp
 - UISHARECollapseOutline
 - UISHAREGetVisible
 - UISHAREQRSCreateProc
 - UISHARESaveAsReadO
 - UISHARETrendBy
 - UITargetCompoundOut
 - UITargetCompoundType

(General) (Declarations)

```
Option Explicit

' Query Report Sheets (QRS) Pivot Table
'
' Purpose: Connect to the QRS database and analyze data
'

' Identifier for particular Query Report Sheet (QRS)

Public QuantitationReportSystemGlobalSetting As String ' global se
Public QRSXMLTable As String ' global se
Public QRSWorksheet As String ' global se
Public QRSXMLMap As String ' global se

Public QRSChartType As XlCh ' Chart Type
Public ChartTitle As ' Chart title

Public PTMWS As Worksheet ' Pivot Table Master Worksheet
Public PTWS As Worksheet ' Pivot Table Worksheet (current)
Public PTCache As PivotCache ' Pivot Table Cache (singular)
Public PTM As PivotTable ' Pivot Table Master (singular)
Public PT As PivotTable ' Pivot Table (current)
Public PTC As Chart ' Pivot Table Chart (current)

Public UIribbon As IRibbonUI ' ribbon
```

Pivot Cache declaration



Microsoft Visual Basic - QRS_TrendCharts_B_04_00.xlsm - [QRSCreatePivot (Code)]

File Edit View Insert Format Debug Run Tools Add-Ins Window Help Type a question for help

Ln 7, Col 34

Project - VBAProject

CreatePivotCache

Pivot Cache creation

```
Sub CreatePivotCache()  
'  
' PivotCache stores summarized information from QRS database  
'  
  
Set PTCache = ActiveWorkbook.PivotCaches.Add(SourceType:=xlDatab  
SourceData:=QRSXMLTab  
  
End Sub
```

Pivot Table creation

```
Sub CreatePivotTable()  
'  
' PivotTable generated without a default setup  
'  
  
Application.DisplayAlerts = False  
  
Set PTMWS = ActiveWorkbook.Worksheets.Add  
  
'Error Condition Possible: You already have a PivotTable and Piv  
PTMWS.Name = "PivotTable" ' !!! IF YOU ARE HERE!!!: Please Delet  
  
Set PTM = PTCache.CreatePivotTable(TableDestination:=PTMWS.Cells  
  
PTM.PreserveFormatting = True  
PTM.SaveData = True
```

Properties - QRSCreatePivot

