

## Investigating the Role of Smad4 In TGF- $\beta$ Signaling using High Density Microarrays

***Collins, P. J.\****, T. Cheung<sup>§</sup>, T. B. Doan\*,  
K. W. Shannon\* and X. Liu<sup>§</sup>\*

Agilent Technologies  
3500 Deer Creek Road, 25U,  
Palo Alto, CA 94304

<sup>§</sup>Department of Chemistry  
and Biochemistry  
UCB215,  
University of Colorado-Boulder  
Boulder, CO 80503

**Synopsis:** Agilent's high-density DNA microarrays are powerful tools in providing answers to complex biological questions. This poster brief presents data from both Human 1 cDNA microarrays and 60-mer oligo microarrays, which illustrate the complexity of the Transforming Growth Factor- $\beta$  (TGF- $\beta$ ) induced transcriptional response in the absence and presence of Smad4. Smad4, a component of the pathway involved in TGF- $\beta$  signal transduction, has been identified as a putative tumor suppressor. Comparison of expression profiles for four tumor cell lines lacking endogenous Smad4 demonstrated that, although considerable heterogeneity exists in their TGF- $\beta$  induced transcriptional responses, different cells share a common set of Smad4-independent TGF- $\beta$  inducible genes. The effect of Smad4 on TGF- $\beta$  signaling was determined by profiling Smad4-transfected colon cancer cells, and both microarray types showed highly correlated data. Northern blotting was used to confirm a number of microarray measurements.




**Agilent Technologies**

Transforming growth factor- $\beta$  is a multifunctional growth factor whose best-known function is to inhibit cell growth and suppress tumor formation. Loss of TGF- $\beta$  growth inhibition is one of the most common cellular events in the pathogenesis of human breast, pancreatic and colon cancers. TGF- $\beta$  signals through a heteromeric signaling complex consisting of Smad2, 3 and 4. Disruption of the Smad signaling complex often leads to tumor formation. We have used both 60-mer oligonucleotide and cDNA microarrays to investigate the role of Smad4 in the TGF- $\beta$  controlled transcription program in tumor cells. These high density DNA microarrays, generated using Agilent's SurePrint inkjet technology, were used to profile global transcriptional regulation in breast, colon and pancreatic Smad4-null tumor cell lines in response to TGF- $\beta$ . Data from both microarray types showed a high degree of correlation in demonstrating that TGF- $\beta$  induces transcriptional activation and repression of genes involved in signal transduction, cell adhesion and transcriptional regulation across the range of cell lines tested. Data from a number of studies are presented comparing expression profiles from Smad4-null tumor cell lines to those from either Smad4-transfected cell lines or normal cell lines. These data indicate that the composition of the Smad signaling complex controls the specificity of TGF- $\beta$  signaling.

- Is Smad4 absolutely required for transmission of TGF- $\beta$  signals?
- Does loss of Smad4 cause similar alterations in TGF- $\beta$  signaling in different tumor cells?
- Is the specificity of TGF- $\beta$  signaling altered in the presence or absence of Smad4?

The Smad4-null tumor cell lines used in this study are listed in **Table 1**. All cell lines were obtained from ATCC and cultured in DMEM supplemented with 10% fetal bovine serum (Invitrogen). Smad4 activity was returned to SW480 cells using retroviral mediated gene transfer. Transfected cells, expressing recombinant Smad4, were selected using FACS sorting as illustrated in Figure **2A**. A western blot was performed using Smad4 antibody after the sorting as shown in Figure **2B**. The response to TGF- $\beta$  stimulation was investigated by treating  $5 \times 10^7$  cells with 100 pM TGF- $\beta$  for 4 hours. After harvesting, total RNA was isolated using RNeasy columns (Qiagen) from untreated and treated cells and labeled using direct incorporation with either cyanine-3 or cyanine-5. Labeled target material representing untreated cells was then mixed with target representing treated cells and the mixture was hybridized to microarrays. Agilent's Human 1 cDNA microarrays (G4100A) were hybridized with cDNA targets which were generated from 10  $\mu$ g of total RNA using Agilent's Fluorescent Direct Label Kit (G2557A). Agilent Human 1A Oligo microarrays (G4110A) were hybridized with cRNA targets which were generated from 50 ng of total RNA using Agilent's Low Input Fluorescent Linear Amplification kit (5184-3523). Microarrays were manufactured using Agilent's SurePrint inkjet technology. Human 1 cDNA microarrays contain more than 12,000 sequence-verified cDNAs from the Incyte Genomics Human UniGene 1 and Human Drug Target clone sets. Human 1A Oligo microarrays consist of experimentally optimized and validated 60-mer oligonucleotide probes, representing more than 17,000 human genes. All experiments were done as dye-swaps with four replicate microarrays being processed for each dye-swap half. Microarrays were scanned using the Agilent dual laser DNA microarray scanner and data extracted using Agilent's Feature Extraction software.



Cell line	Origin
MDA-MB468	Breast cancer cells
SW480	Colon cancer cells
CFPAC-1	Pancreatic cancer cells
Hs 766 T	Pancreatic cancer cells



**Figure 2A.** Selection of cells expressing recombinant Smad4




**Figure 2B.** Western Blot image using Smad4 antibody




**Figure 3: Expression profiles of Smad4-null tumors** of SW480 (A), Hs 766 T (B), MDA-MB468 (C) (polarity). Data points represent Resolver-comb against signal intensities. The Agilent Feature each experiment. Log ratios colored blue are u (significantly greater than 0) and those in green

Figure 3E presents a comparison of the number of genes different from 0 for SW480 and CFPAC-1 than for the other cell lines. Genes were identified for all cell lines which had




**Figure 4: Measuring the effect of Smad4 on e**  
The effect of Smad4 on global gene expression  
both Human 1 cDNA microarrays (**A**) and Hum  
significantly up-regulated in Smad4-positive ce  
Although the resulting profiles show a low leve  
by both microarray platforms. The differentia  
regulated in the Smad4 expressing cells relativ  
gene presented, Inhibitor of DNA binding 1 (**D**)




**Tumor cell lines treated with TGF- $\beta$ .** The effect on gene expression of 4 h treatment with TGF- $\beta$  in each cell line (A, B, C, and D) was measured using 8 replicate Agilent Human 1 cDNA microarrays (4 in each channel) and 8 replicate Agilent Human 1A Oligo microarrays (4 in each channel). The extraction error model was utilized by Resolver to determine log ratio value significances ( $p < 0.001$ ) for differentially expressed gene elements (not significantly different than 0), those shown in red are up-regulated and those shown in green are down-regulated (significantly less than 0).

For the identification of TGF- $\beta$  responsive genes in the four cell lines. More log ratios were determined to be significantly different from 0 for Hs7667 T and MDA-MB468, indicating overall differences in responsiveness to TGF- $\beta$ . However, these genes have previously been reported to be regulated by TGF- $\beta$ .



**Expression profiles in the absence of TGF- $\beta$  using Human 1 cDNA and Human 1A Oligo microarrays.** The effect on gene expression of 4 h treatment with TGF- $\beta$  in each cell line (A, B, C, and D) was measured using 8 replicate Agilent Human 1 cDNA microarrays (4 in each channel) and 8 replicate Agilent Human 1A Oligo microarrays (4 in each channel). The extraction error model was utilized by Resolver to determine log ratio value significances ( $p < 0.001$ ) for differentially expressed gene elements (not significantly different than 0), those shown in red are up-regulated and those shown in green are down-regulated (significantly less than 0). For the identification of TGF- $\beta$  responsive genes in the four cell lines. More log ratios were determined to be significantly different from 0 for Hs7667 T and MDA-MB468, indicating overall differences in responsiveness to TGF- $\beta$ . However, these genes have previously been reported to be regulated by TGF- $\beta$ .

**Figure 5: Identification of Smad4- dependent genes in colon cancer cells using both Human 1 cDNA microarrays and Human 1A Oligo microarrays.** The response of Smad4-transfected cells to TGF- $\beta$  stimulation was determined by hybridizing TGF- $\beta$  treated cells against untreated cells to both Human 1 cDNA microarrays (A) and to Human 1A microarrays (B). Expression data from both microarray platforms correlate very highly, with less than 0.4% of genes being anticorrelated. Differential expression measurements from both microarray types corresponding to a number of previously well-characterized TGF- $\beta$  responsive genes are listed in Table 2. Data from this study indicate that, as well as being TGF- $\beta$  independent, expression of these genes is Smad4- dependent.



**Table 2. Examples of TGF- $\beta$  -Dependent Smad4-Dependent Genes**

Gene	Response	
	Human 1 cDNA	60-mer Oligonucleotide
junB proto-oncogene	2.3-fold Up	2.7-fold Up
Smad7	2.5-fold Up	1.8-fold Up
Lamin alpha 3	3.2-fold Up	4.9-fold Up
Disintegrin/Zinc Metalloprotease	2.0-fold Down	1.6-fold Down

### Conclusions

- Smad4 is not absolutely required for TGF- $\beta$  induced signaling but does alter the specificity of the response.
- Different tumor cells exhibit considerable heterogeneity in their TGF- $\beta$  induced transcriptional response but share a common set of Smad4-independent TGF- $\beta$  inducible genes.
- The precision of Agilent's ink-jet microarrays enables accurate characterization of the TGF- $\beta$  induced transcription response.



## Ordering Information

[www.agilent.com/chem/dna](http://www.agilent.com/chem/dna)  
u.s. and canada 1 800 227 9770  
japan +0120 477 111  
europe: [marcom\\_center@agilent.com](mailto:marcom_center@agilent.com)  
global: [dna\\_microarray@agilent.com](mailto:dna_microarray@agilent.com)

© Agilent Technologies, Inc. 2003

Research Use Only

LifeSeq is a trademark or registered trademark of Incyte Genomics, Inc. in the U.S. and other countries. Rosetta Resolver is a U.S. registered trademark of Rosetta Inpharmatics. Information, descriptions and specifications in this publication are subject to change without notice. Printed in the U.S.A.

Printed in the U.S.A.  
May 1, 2003  
5988-9292EN



**Agilent Technologies**