

X-RAY SCINTILLATOR

ACS ALS FOS




ACS[®] Amorphous-Carbon Plate
with CsI Scintillator

ALS[®] Aluminum Plate
with CsI Scintillator

FOS[®] Fiber Optic Plate
with CsI Scintillator

Scintillator: CsI(Tl)
Columnar Structure





HAMAMATSU

X-RAY SCINTILLATOR

ACS | ALS | FOS

STRUCTURE



FEATURES

- **Large Format** Maximum size: 468 mm (17") x 468 mm (17") for ACS and ALS
- **High Light Output** 2.5 times higher with ACS-HL type (CsI 600 μm) than Lanex-R (powdery phosphor). [Typ.]
- **High Resolution** 20 Lp/mm at CTF 13 % FOS-HR type (CsI 150 μm). [Typ.]


SELECTION GUIDE

Product name	Structure	Availability in dimension			Scintillator Thickness (μm)	Features	Applications
		Scintillator effective area (mm)		Substrate Thickness (mm)			
		Max.	Min.				
ACS	Amorphous-Carbon Plate with CsI Scintillator	440 × 440 (17" × 17")	14 × 14	0.5 or 2	600 Max.	High resolution, Large format	Dental inter oral, Mammography, Chest examination
ALS	Aluminum Plate with CsI Scintillator	440 × 440 (17" × 17")	14 × 14	1		High light output, Large format	Dental-panoramic, Chest examination
FOS	Fiber Optic Plate with CsI Scintillator	240 × 180 (9" × 7")	10 × 10	1 to 3		X-ray shield, Low energy X-ray detection	Dental inter oral, Dental-panoramic, Mammography


APPLICATION EXAMPLE



Dental



Mammography



Chest examination

STANDARD PRODUCTS LINE-UP


* Available in various sizes and scintillator thickness.

	Type No.	Scintillator type [Ⓐ]	Outer dimension (mm)	Effective area (mm)	Substrate Thickness (mm)	CsI Thickness (μm)	Relative light output [Ⓑ] (% Typ.)	CTF [Ⓒ] (% Typ.)	Type [Ⓓ]
FOS	J6671	CsI (TI)	30.5 × 21	27 × 17	3	150	70	18 [Ⓔ]	HL
	J6671-01					150	40	33 [Ⓔ]	HR
	J6673	CsI (TI)	50 × 10	47 × 7	3	150	70	18 [Ⓔ]	HL
	J6673-01					150	40	33 [Ⓔ]	HR
	J6675	CsI (TI)	18 × 18	15 × 15	3	150	70	18 [Ⓔ]	HL
	J6675-01					150	40	33 [Ⓔ]	HR
	J6677	CsI (TI)	50 × 50	47 × 47	3	150	70	18 [Ⓔ]	HL
	J6677-01					150	40	33 [Ⓔ]	HR
J6679	CsI (TI)	φ26.5	φ23.5	3	150	70	18 [Ⓔ]	HL	
J6679-01					150	40	33 [Ⓔ]	HR	
ACS	J8734	CsI (TI)	50 × 50	48 × 48	0.5	150	125	12 [Ⓔ]	HL
	J8734-01					150	50	25 [Ⓔ]	HR
	J8977	CsI (TI)	468 × 468	440 × 440	2	600	250	33 [Ⓕ]	HL
ALS	J8978	CsI (TI)	50 × 50	48 × 48	1	150	70	20 [Ⓔ]	
	J9857	CsI (TI)	468 × 468	440 × 440	1	600	150	35 [Ⓕ]	



- [Ⓐ]150 μm thickness for the CsI(Tl) scintillator.
- [Ⓑ]Relative values, with 100 % being equal to the light output from conventional phosphor screen (Lanex-R). Light output was measured by CCD with lens coupling under the following conditions : (X-ray tube voltage 60 kV p, aluminum filter 1 mm thick)
- [Ⓒ]CTF (contrast transfer function) CsI(Tl) : X-ray tube voltage 60 kV p, aluminum filter 1 mm thick
- [Ⓓ]HL: high light output type, HR: high resolution type
- [Ⓔ]at 10 lp/mm [Ⓕ]at 3 lp/mm

X-RAY TRANSMITTANCE



* a-C: Amorphous-Carbon
Amorphous Carbon has good X-ray transmittance characteristics because it is a light element material. In addition, it is a glass like material with no particle causing blemish defects. It can be polished to a good of flatness for combination with a 2D image device. Also, it is environmentally safe.

WHAT'S FOP?


Fiber Optic Plate: FOP

The FOP is an optical device consisting of millions of glass fibers of several micrometers in diameter, bundled parallel to one another. Since light is transmitted through each fiber, an image appears to float. The image can be transferred from one end of the fiber to the other without any distortion. FOPs are widely used as optical devices that replace optical lens.



The reason is the "optical fiber structure".

Why does the image appear to float?



X-RAY RELATED PRODUCTS

● MICROFOCUS X-RAY SOURCE

Due to the minute size of the focal spot, Hamamatsu microfocus X-ray sources allow capturing clear, sharp images even if magnified. Various types of microfocus X-ray sources are available including the sealed-off types of 80 kV, 90 kV, 100 kV, 110 kV, 130 kV and 150 kV and the open types of 160 kV.



● X-CUBE™ (COMPACT X-RAY CCD CAMERA)

X-CUBE™ is a compact X-ray CCD camera designed for non-destructive inspection. Using a general-purpose CCD chip mounted in a rugged but lightweight camera head, X-CUBE™ makes X-ray imaging as easy as handling ordinary CCD camera.



● X-RAY SHIELD FIBER OPTIC PLATE

XRS-FOPs provide enhanced X-ray absorption rate about 50 times higher than normal FOPs. Almost all X-rays which have penetrated the scintillator and have not been converted into light are absorbed in the XRS-FOP. This eliminates X-ray damage of image sensors such as CCDs. XRS-FOPs are available in sizes up to 50 mm × 50 mm



Type No.	Fiber diameter	Numerical aperture (N.A.)	Resolution (Lp/mm)	Absorption material	Thermal expansion coefficient ($\times 10^{-7}/^{\circ}\text{C}$)	Transmittance (%)	Dimensions (mm)	Thickness (mm)
J6745-01	6 μm	1.0	90.5	Included	96	60	30 × 21	3

*PATENT: USA: USP6531225, USP6762420, Europe: EP1024374B, Japan: 1832818, 3126715, 3566926, China: ZL99801885.6

Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult with our sales office. Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2009 Hamamatsu Photonics K.K.

HAMAMATSU

WEB SITE www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, P. O. Box 6910, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road Welwyn Garden City Hertfordshire AL7 1BW, United Kingdom, Telephone: 44-(0)1707-294888, Fax: 44(0)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171-41 SOLNA, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia: S.R.L.: Strada della Moia, 1/E, 20020 Arese, (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741 E-mail: info@hamamatsu.it

TMCP1031E04
JUN. 2009 IP
Printed in Japan (1000)