

# External Quantum Efficiency Measurement System C9920-12

External quantum efficiency measurement system for light emitting device utilizing an integrating sphere



## Luminous Efficiency can be measured precisely without Effects of Emission Angle Distribution Characteristics.

The C9920-12 External Quantum Efficiency Measurement System is a system for measuring the external quantum efficiency of a light emitting device such as organic/inorganic LED and LED by electro-luminescence (EL) method.

The luminous efficiency of the light emitting devices is related to various elements, such as the absorption of the light emitting part and glass substrate, and the reflectance of the reflective mirror surface, so on. The C9920-12 can measure the efficiency of external emissions for the applied current (voltage), inclusive of these elements related to efficiency. The current (voltage) applied to the sample can be set on



the software. It is applied to the sample in steps from the minimum value to the maximum value for the current (voltage) set in advance from a power source meter. With an integrating sphere, highly precise emission efficiency measurements can be made independently of the emission angle distribution characteristics.

Employing Hamamatsu multi-channel detector, the spectrum can be measured instantaneously at each step. The spectral data from each step is recorded by the software, and the external quantum efficiency, I-V-L characteristics and chromaticity can be displayed.



## External quantum efficiency measurement system utilizing an integrating sphere.

The luminous efficiency of the light emitting device including the optical effect from the glass substrate and reflective mirror, etc. is measured rather than just that of the emitting part. The power source meter is controlled on software and the spectrum at each current (voltage) step is measured automatically.



### FEATURES

- With an integrating sphere the measurement does not influenced on the emission angle characteristics of the sample
- The power source meter (KEITHLEY 2400) can be controlled on the software
- The spectrum can be measured instantaneously with Hamamatsu multi-channel detector for each step of the applied current (voltage)
- Using a high-performance, BT(buck thinned) cooled CCD ,the system has ultra-high sensitivity
- Settings, measurements and displays of results are simplified by the dedicated software
- On the software, it is possible to display various types of graphs, such as current/current density/voltage-spectrum, current/current density/luminous efficiency or current/current density/voltage-chromaticity.
- The system can be easily extended to absolute PL quantum yield measurements and brightness and light distribution characteristics measurement systems with additional parts

### SPECIFICATIONS


Integrating sphere	3.3 inch inside diameter, reflective material: Spectralon
Detector	BT(buck thinned) cooled CCD
Cooling temperature (CCD)	-15 °C
No. of photosensitive device channels	1024 ch
Wavelength measurement range	380 nm to 780 nm (detector: 200 nm to 950 nm)
Fiber probe	1.5 m
Range of measurement for luminous flux	0.00013 lm to 0.12 lm (with an emission area 2 mm square and white light)

# Simple and easy operativity is realized with the dedicated analysis software.

Software for measuring external quantum efficiency

Measurements for external quantum efficiency of light emitting device can be made on the software.

Screen displaying measurement results




Minimum and maximum value of current/voltage value, and current/voltage value steps and the like can be set in a window at the upper left of the display screen. Raw measurement data (emission spectrum etc.) is displayed in the window. Sample monitoring measurements, data during a measurement and the like are displayed. Various graphs can be displayed in the figures on the right side of the display screen.


Graph display

The results of measurements can be displayed in various graph formats.


Current density/Voltage




Current density/Chromaticity coordinates




Current density/Dominant wavelength




Current density/Excitation purity



Current density/Power efficiency



Wavelength/Intensity



Current density/External quantum efficiency

Time/emission spectrum

Time/Current (voltage)

Time/Chromaticity coordinates

Current/external differential quantum efficiency

Current density/Luminous efficiency (cd/A)

\* Voltage can also be displayed in place of current density.

# A wide variety of measurement targets can be accommodated by adding optional parts and components.

## OPTIONS



### 150 W CW Xenon Light Source with Input Optics L10092

This is necessary when supplementing for sample absorption. It can be used jointly with an organic LED quantum absorption measurement device.



### Integrating Sphere Adaptor A10093

This is an adaptor for connecting an absorption supplementing light source to the integrating sphere.




### EL Measurement Sample Holder A9924-04

This is a sample holder that uses pin type connectors.



### EL Measurement Sample Holder A9924-06

This is a general purpose sample holder that uses alligator clips. A light shielding cover is included.



## Related Products

There are many products for combined configurations in the C9920 series, and expansion to the following products is easy.

### Absolute PL Quantum Yield Measurement System C9920-01,-02

The photo-luminescence quantum yield is measured for photo-luminescent materials using photoluminescence (PL) method.


Excitation at various wavelengths is possible with a band-pass filter or a monochromator. Various sample holders are provided, so measurements can be made not only on thin films, but also on solutions and powders.

### Brightness Light Distribution Characteristics Measurement System C9920-11


Measurements such as brightness for each emission angle, emission spectrum and chromaticity are possible using a rotating stage.

The light emitting device emission brightness, spectrum and emission angle distribution are measured for each of the angle steps that has been set.

## Block Diagram



## Block Diagram



- ★ DynaSpect is a registered trademark of Hamamatsu Photonics K.K.
  - Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult with our sales office.
  - Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.
- © 2006 Hamamatsu Photonics K.K.

# HAMAMATSU

[www.hamamatsu.com](http://www.hamamatsu.com)

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: [export@sys.hpk.co.jp](mailto:export@sys.hpk.co.jp)

U.S.A. and Canada: Hamamatsu Corporation, Systems Division: 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1) 908-231-0960, Fax: (1) 908-231-0852, E-mail: [usa@hamamatsu.com](mailto:usa@hamamatsu.com)

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658, E-mail: [info@hamamatsu.de](mailto:info@hamamatsu.de)

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10, E-mail: [infos@hamamatsu.fr](mailto:infos@hamamatsu.fr)

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, U.K., Telephone: (44) 1707-294888, Fax: (44) 1707-325777, E-mail: [info@hamamatsu.co.uk](mailto:info@hamamatsu.co.uk)

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171-41 Solna, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01, E-mail: [info@hamamatsu.se](mailto:info@hamamatsu.se)

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1/E 20020 Arese (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741, E-mail: [info@hamamatsu.it](mailto:info@hamamatsu.it)

Cat. No. SSMS0018E03  
APR/2008 (2006.10) HPK  
Created in Japan