

DynaSpect® Series

DynaSpect is the name for a series of spectrophotometer devices.

Photonic Multichannel Analyzer

PMA-12

Scientific applications

UV to visible spectroscopy

Fluorescence spectroscopy

Raman scattering

Chemiluminescence analysis

Liquid chromatography

Gas chromatography

ICP emission analysis

Discharge spectrum analysis

Combustion analysis

Micro spectroscopy

Industrial applications

Water quality testing

Evaluation of light emitting devices and light sources

Chromaticity measurements

Impurities testing

Film thickness measurements

UV radiation measurements

Plasma monitors

Fruit testers

Combustion monitoring

Color filter evaluation

HAMAMATSU

PMA Photonic Multichannel Analyzer

**A compact unit containing a spectrometer, photo-detector and power supply.
Use of an optical fiber input makes spectral measurements easy.**

It becomes easier to use of PMA-12 because of USB2.0 interface



The PMA-12 is compact spectral measurement apparatus that combines a spectrometer and optical detector into one unit. An optical fiber is used. Because of the high sensitivity, spectra can be obtained easily just by bringing the optical fiber close to the sample in normal applications, without a special light collection system. Since the spectrometer and photo-detector are fixed, the PMA-12 is stable and can be used with confidence for long periods of time. The wavelength axis and spectral response characteristics are already calibrated, so spectral measurements can be carried out easily and accurately.

C10025-01 High sensitivity superior cost-performance model
C10544-01 C10544-02

C10025-01, C10544-01, -02, which have the thermoelectric cooling type CCD linear image sensor that is used for astronomical observation, have realized both high performance and low price by rational design. The wavelength range for measurement is 300 nm to 800 nm for the C10025-01/C10544-01 and 340 nm to 830 nm for the C10544-02.

C10028-01 C10028-02 Near infrared model

These are models using InGaAs linear image sensors and capable of measurements of reflection and absorption spectra in the near infrared with a large dynamic range. The wavelength range for measurements is 900 nm to 1650 nm for the C10028-01 and 1600 nm to 2350 nm for the C10028-02.

C10027-01 C10027-02 Ultra-high sensitivity model



This model uses a thermoelectrically cooled, back-thinned CCD linear image sensor that with higher sensitivity and lower noise than the C10025-01. The C10027-01 is an ultra-high sensitivity model that combines this sensor with a small Czerny-Turner spectrograph capable of measurements over a wide range from the ultraviolet to the near infrared with high wavelength resolution. The wavelength range for measurements is 200 nm to 950 nm for the C10027-01 and 350 nm to 1100 nm for the C10027-02.

C10029-01 High time resolution model

Coupling an image intensifier with a thermoelectrically cooled, back-thinned CCD linear image sensor, it is possible to have both high-speed gate measurements at a maximum of 10 ns and ultra-high sensitivity. This model is capable of high temporal resolution measurements in the nanosecond range and measurements of faint light.


FEATURES

- Spectrometer, photo-detector and power supply in a compact unit
- Real-time measurements (Simultaneous measurement of multiple wavelengths possible)
- Easy measurements with optical fiber
- Spectral response and wavelength axis characteristics calibrated
- Wide range of variations



MEASUREMENT MODES

- Standard measurements**
 This is the most basic measurement mode.
 Applications: emission spectra for light sources, fluorescence, plasma and the like.
- Reflective measurements**
 This is the measurement mode for finding spectral reflectance.
 Applications: reflectance measurements for optical filters, coatings and the like.
- Transmittance and absorption measurements**
 This is the measurement mode for finding spectral transmittance and absorption.
 Applications: measurements of transmittance and absorption in optical filters, films, solutions and the like.
- Chromaticity measurements (light-source color)**
 This is the measurement mode for finding the light-source color for luminous bodies.
 Applications: color evaluation in light sources for illumination, LEDs and the like.
- Chromaticity measurements (object color)**
 This is the mode for finding the color of objects that are either reflective or transmit light.
 Applications: color evaluation of paint, fabric, printed matter and the like.



DISPLAY MODES

Spectrum display

Counts Intensity

Wavelength nm

Display of changes over time

Counts Intensity

Time

3-D display

Wavelength nm

Time

Reflectivity display

% Percent

Wavelength nm

Transmittance display

% Percent

Wavelength nm

Absorbance display (OD)

OD

Wavelength nm

Color coordinate display

Analyze [XY]

Angle(Degree) 2 10

Information

Tristimulus Values

X: 192109

Y: 300257

Z: 200651

Chromaticity Coordinates

x: 0.323944

y: 0.641209

Dominant Wavelength: 552.99 nm

Excitation Purity: 90.792 %

Spatial color coordinate display

Analyze [Lab]

Angle(Degree) 2 10

Information

Tristimulus Values

L: 144134

a: 160452

b: 445023

CIE(LAB Color Space)

L: 13.2971


a: -39.8034

b: -16.943

SYSTEM EXAMPLES

Light source measurements

Measurement of emission spectra in light sources such as lamps and LEDs




<Configuration>

- Standard PMA-12 configuration (C10025, C10544, C10027, etc.)
- C10471 PMA data analyzer


<Applications>

- Evaluation of color temperature and color rendering properties in light sources for illumination
- LED chromaticity evaluations
- Special applications of light source spectral evaluations


Germicidal lamp emission spectrum



Metal halide lamp emission spectrum




White LED emission spectrum




■ Analysis of light source color by emission spectrum (chromaticity, color temperature, color rendering properties, etc., possible)

Metal halide lamp chromaticity evaluation




White LED chromaticity evaluations



Emission spectrum measurements


Emission spectrum measurements for plasma, electric discharge, ablation and the like




<Configuration>

- Standard PMA-12 configuration (C10027, C10029, etc.)
- C10471 PMA data analyzer
- A8482 microscope position checking optics
- A4869 UV lens
- DG535 delay generator (for gate operation)


Electric discharge emission spectrum



Emission spectrum during oxide film etching




Emission spectrum during nitride film etching



<Applications>



- Plasma component analysis
- Analysis of various emission phenomena

Temporal changes in plasma emission spectra




Example of temporal resolution (gate operation) measurements

PMA-12 (C10029-01)




Reflective spectrum measurements

Measurement of spectral reflectance in optical filters, anti-reflective films (AR coatings) and the like




- <Configuration>**
- Standard PMA-12 configuration (C10025, C10544, C10027, etc.)
 - C10471 PMA data analyzer
 - L6759 high stability 150 W xenon light source
- <Applications>**
- Inspection of coatings
 - Monitoring thin film growth

AR coating reflection spectrum




Object color measurements

Object color evaluation of paint, fabric, printed matter and the like




- <Configuration>**
- Standard PMA-12 configuration (C10025, C10544, C10027, etc.)
 - C10471 PMA data analyzer
 - L6758 halogen light source
- <Applications>**
- Paint inspections
 - Color evaluations in printed matter, fabric, plastics, etc.

Paper object color (chromaticity coordinates)




Absorption spectrum measurements

Spectral transmittance and absorption measurements in optical filters, films, solutions and the like




- <Configuration>**
- Standard PMA-12 configuration (C10025, C10544, C10027, etc.)
 - C10471 PMA data analyzer
 - L6759 high stability 150 W xenon light source
 - A6751 solution sample holder for transmissivity and fluorescence measurements
- <Applications>**
- Absorption spectrum evaluations for solutions and films
 - Component analysis for samples
 - Monitoring chemical changes


Component analysis of plastics using transmission spectra (polycarbonate and PET resins)




MMA and PMMA transmission spectra



Didymium film absorption spectrum




Changes in transmissivity in the polymerization from MMA to PMMA (wavelength: 1615 nm)



Microscopic spectral measurements

Spectral distribution measurements under a microscope




- <Configuration>**
- Standard PMA-12 configuration (C10027, C10029, etc.)
 - C10471 PMA data analyzer
 - A8482 C-mount position checking optics
- <Applications>**
- Measurement of bioluminescence
 - Measurements on semiconductor wafer, LCD and other microstructures


SYSTEM EXAMPLES

Emission spectrum measurements


For fluorescent samples such as fluorescent lamps and EL devices




Fluorescence indicator (Fluorescein) emission spectrum



Chemiluminescence emission spectrum



Fluorescent lamp fluorescent body emission spectrum



<Configuration>


- Standard PMA-12 configuration (C10025, C10544, C10027, etc.)
- C10471 PMA data analyzer
- Excitation light source: laser, xenon lamp, etc.
- A6751 solution sample holder for transmissivity and fluorescence measurements

<Applications>

- Fluorescence spectroscopy
- Monitoring chemical light emissions

Film thickness measurements


Film thickness measurements using spectral reflectance or transmittance



<Applications>

- Monitoring thin film growth
- Film thickness management
- Resist film thickness measurements

ITO film interference spectrum



<Configuration>


- Standard PMA-12 configuration (C10025, C10544, C10027, etc.)
- C10471 PMA data analyzer
- L6759 high stability 150 W xenon light source

Optical NanoGauge (Interferometric Film Thickness Measurement System) C10178, C10323


We can offer specialized machine for film thickness measurement. Please refer the detail in specific brochure.

Quantum yield measurement system

Measurement of quantum yield, external quantum efficiency, brightness light distribution characteristics



[Screen showing emission spectrum]



<Applications>

- Research of fluorescence materials in physics or chemistry
- Quantum yield measurement of emission materials
- Internal quantum yield measurement of fluorescence materials

Absolute PL Quantum Yield Measurement System C9920-01, -02

External Quantum Efficiency Measurement System C9920-12

Brightness Light Distribution Characteristic Measurement System C9920-11

We can offer specialized machine for OLED measurement. Please refer the detail in specific brochure.

OPTIONS



Solution sample holder for transmissivity and fluorescence measurements A6751

This is a dedicated holder with an integrated condensing lens for the use with vials.



Reflection measurement optics A9665

These are optics making it possible to illuminate the sample at 45° from the light source and measure the reflected light.



Variable angle reflection measurement optics A10687

These are optics making it possible to change the angle of input and output ports at maximum 65° and measure the reflected light and fluorescence.



Digital delay generator DG535

This outputs the gate pulse necessary for an external trigger and gate operation.



2 split fiber A10193-01,-02

It is very useful for reflectance measurement or film thickness measurement. We have two kind of fiber. One is A10193-01 for visible range and the other is A10193-02 for from visible range to near infrared range.



C-mount fiber adapter A6399

This is an adapter for securing the fiber input optics to the C-mount of a microscope or the like.



C-mount position checking optics A8482

In addition to the function of the C-mount fiber adapter, the measurement position can be checked.



C-mount position checking fiber adapter A9607

In addition to the function of the C-mount fiber adapter, the measurement position can be checked. Measurements in the UV range are possible.



PMA data analyzer C10471

A data analyzer is provided. There are the C10471-01 notebook model and the C10471-02 desktop model.



UV lens A4869

Condensing lens for UV. f=50 mm, F3.5 (A6399, A8482 required)



Integrating sphere A5640

This is the integrating sphere for getting complete diffuse light. You can get even intensity light without spread of light source or influence of directional characteristics.



Fading light adaptor A10474-01

This adaptor is used when the light power is too strong. It can reduce the input light power by using a pin-hole.

(fading rate approx 1/20 to 1/500)



Light source with monochromator

This light source can choose the excitation wavelength by using the monochromator.

(selectable wavelength area 250 nm to 800 nm)



Halogen light source L6758

This is a halogen light source with output wavelengths from 400 nm to 1000 nm for excitation and absorption measurements.

* Light guide connector A10194-01 is needed to connect with 2 split fiber.



High stability 150 W xenon light source L6759

This is a high stability xenon light source with output wavelengths from 250 nm to 1000 nm for excitation and absorption measurements.



Light source with wide wavelength range L8128-01




This is a light source with output wavelengths from 200 nm to 950 nm and no bright lines for transmissivity, absorption and reflection measurements.

PMA software library U10472-01

This is the software library which controls the PMA-12 series.

PMA color measurement library U10473-01

This is the software library which controls the PMA-12 series and calculates the chromaticity.



SPECIFICATIONS


Model	C10025-01	C10544-01	C10544-02	C10027-01	C10027-02	C10028-01	C10028-02	C10029-01
Photo-detector	CCD linear image sensor			BT- CCD linear image sensor		InGaAs linear image sensor		I.I. + BT- CCD linear image sensor
Wavelength (nm)	300 to 800		340 to 830	200 to 950	350 to 1100	900 to 1650	1600 to 2350	200 to 860
Wavelength resolution (FWHM) ¹	< 3 nm			< 2 nm	< 2.5 nm	< 9 nm		< 3 nm
Exposure time	18 ms to 64 s	19 ms to 64 s		19 ms to 64 s		5 ms to 64 s	5 ms to 50 ms	19 ms to 64 s
Gate time ²	-			-		-	-	≥ 10 ns
Gate repetition	-			-		-	-	≤ 200 kHz
Number of photosensitized device channels	1024 ch			1024 ch		256 ch		900 ch
Pixel size	24 μm × 1.53 mm	24 μm × 3.07 mm		24 μm × 2.928 mm		50 μm × 250 μm		24 μm × 2.928 mm ³
Device cooling temperature	0 °C			-15 °C		-10 °C		-15 °C ³
Read-out noise (electrons)	18			16		12 500		16 ³
Dark current (electrons/scan)	200 (0 °C : 20 ms)	400 (0 °C : 20 ms)		75 (-15 °C : 20 ms)		20 000 (-10 °C : 20 ms)	2.5 × 10 ⁷ (-10 °C : 20 ms)	75 ³ (-15 °C : 20 ms)
AD resolution	16 bit							
Spectrograph	Concave spherical grating type			Czerny-Turner type				
Spectrograph F number	3			4				
Fiber receiving area	φ800 μm			φ1 mm				
Fiber type	Single fiber · SUS tube			Bundled fiber φ12 mm SUS tube				
Fiber length	1.5 m							
External trigger input	TTL level / High impedance							
Interface	USB 2.0							
Power supply	AC 100 V to AC 240 V, 50 Hz / 60 Hz (Power supply voltage variation ±10 %)							

¹ Confirmed with mercury and argon atomic beams.
² The gate time is controlled by the external gate pulse width.
³ I.I. characteristics are not included.


DIMENSIONAL OUTLINES (Unit :mm)

● Main unit


• C10025-01, C10544-01, C10544-02




• C10027-01, -02



• C10028-01, -02




• C10029-01




● Fiber input optics for C10025

(approx. 86 g)



● Fiber input optics for C10544, C10027, C10028, C10029

(approx. 100 g)



PMA SOFTWARE U6039-01

- Measurement functions Monitoring measurement
Data measurement
- Temporal resolution measurement functions ... Temporal fluctuation of spectra
Temporal fluctuation in reflectivity and transmissivity
- Data acquisition condition settings Exposure time settings
Memory integration count assignment
- Calibration/correction Wavelength axis calibration
Sensitivity inconsistency calibration
Dark current correction
- Display functions Spectrum display
Display temporal waveform fluctuations
- Wavelength axis display Wavelength, wave number, Raman shift, energy (eV)
- Brightness axis display Linear, logarithm
- Cursor analysis functions Wavelength (wave number, etc.) vs. intensity
Peak detection
FWHM measurement
Integrated intensity
- Other analytical functions Smoothing
Differential waveform
Color calculation (XYZ, xy, uv, Lab)

★ DynaSpect is a registered trademark of Hamamatsu Photonics K.K.
 ★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.
 ★ The spectral response graphs in this document are representative examples, so they are not guaranteed.
 • Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.
 Specifications and external appearance are subject to change without notice.
 © 2008 HAMAMATSU PHOTONICS K.K.

HAMAMATSU

Homepage Address www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Systems Division
 812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A. and Canada: Hamamatsu Corporation, Systems Division: 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1) 908-231-0960, Fax: (1) 908-231-0852, E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658, E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10, E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road Welwyn Garden City Hertfordshire AL7 1BW U.K., Telephone: (44)1707-294888, Fax: (44)1707-325777, E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171-41 Solna, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01, E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1/E 20020 Arese (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741, E-mail: info@hamamatsu.it

Cat. No. SDSS 0008E05
 APR/2008 (2008.2) HPK
 Printed in Japan