

## PSV 9.1 Scanning Vibrometer Software

Polytec's Scanning Vibrometer Software is the "brain" behind all PSV systems. The software guides and controls all measurement tasks. This includes setting measurement parameters and scanning the laser as well as acquiring, displaying, processing and exporting measurement data. The latest release, PSV 9.1, incorporates the most recent recommendations of the vibrometer user community, leading to improved system functionality and performance.

### How Can You Benefit from the Latest Innovations?

Polytec will design your personal software maintenance plan entitling you to receive your copy of the latest release.

### Further Information

Please consult the documentation delivered with your software, including release notes, for a detailed description of

- further improvements
- hardware requirements and compatibility
- supported software interfaces and bug fixes



### New Features

- Automatic 3D Alignment
- Improved High Contrast Laser Display
- Extended UNDO / REDO in Acquisition Mode
- TEDS support
- Calculation of Damping Values
- Differential cursor: frequency readout
- 64-bit support
- Compensation in Time Domain

# PSV 9.1 Scanning Vibrometer Software

## Setup, Measurement and Data Analysis

### Datasheet



# New Features in PSV 9.1

## **Automatic 3D Alignment**

Machine vision is used to maneuver the 2 remaining lasers to one reference point indicated by 1 laser. This procedure works for all set-ups on smooth surfaces. Switching the lasers is no longer necessary. Only the TOP laser is moved manually.



### **Your Benefit:**

- Speeds up the alignment process by a factor of 3

## **Improved High Contrast Laser Display**

The High-Contrast Laser Display improves the visibility and clarity of the laser beam on the live video image. The laser position is now tracked and graphically marked, improving the visibility and operability.



### **Your Benefit:**

- Improved alignment precision and point definition with shiny objects.

## **Extended UNDO / REDO in Acquisition Mode**

UNDO is now available for acquisition settings, alignments, geometry and camera settings.



### **Your Benefit:**

- Accelerates the setup process

## **TEDS Support**

Support of TEDS sensors (Transducer Electronic Data Sheet; IEEE 1451). The sensitivity specification of each reference sensor is automatically stored in the A/D settings when connected to a PSV-F-500 Front-End.



### **Your Benefit:**

- Removes a source of error
- Accelerates the setup process

### Calculating the Damping Factor

For damping estimations the Q factor and Zeta are added.



#### Your Benefit:

- Easy estimation of modal damping parameters for an initial validation of FE results

### Differential Cursor: Frequency Readout

Frequencies are calculated from two cursors that can be positioned anywhere on a time trace. Waveform frequency content can be quickly and easily studied.



#### Your Benefit:

- Quicker evaluation in time domain

### Compensation in Time Domain

The propagation delay of the decoder relative to a reference channel is now compensated in time domain. In the past the compensation was only applied in the frequency domain.



#### Your Benefit:

- Timing between laser vibrometer and reference sensor can be used for analysis

### 64-bit Support

PSV software is running as a 64-bit application in a 64-bit environment. This allows access to data from other 64-bit applications via Polytec File Access (PFA). Over 4 GBytes of RAM is used.



#### Your Benefit:

- Enhances speed for large data files
- Improved overall performance in acquisition and presentation modes



 **Polytec GmbH  
(Germany)**  
Polytec-Platz 1-7  
76337 Waldbronn  
Tel. +49 7243 604-0  
info@polytec.de

**Polytec GmbH  
(Germany)**  
**Vertriebs- und  
Beratungsbüro**  
Schwarzschildstraße 1  
12489 Berlin  
Tel. +49 30 6392-5140

 **Polytec, Inc.  
(USA)**  
North American  
Headquarters  
16400 Bake Parkway  
Suites 150 & 200  
Irvine, CA 92618  
Tel. +1 949 943-3033  
info@polytec.com

**Central Office**  
1046 Baker Road  
Dexter, MI 48130  
Tel. +1 734 253-9428

**East Coast Office**  
25 South Street, Suite A  
Hopkinton, MA 01748  
Tel. +1 508 417-1040

 **Polytec Ltd.  
(Great Britain)**  
Lambda House  
Batford Mill  
Harpenden, Herts AL5 5BZ  
Tel. +44 1582 711670  
info@polytec-ltd.co.uk

 **Polytec France S.A.S.**  
Bâtiment Orion – 1er étage  
39, rue Louveau  
92320 Châtillon  
Tel. +33 1 496569-00  
info@polytec.fr

 **Polytec Japan**  
Arena Tower, 13th floor  
3-1-9, Shinyokohama  
Kohoku-ku, Yokohama-shi  
Kanagawa 222-0033  
Tel. +81 45 478-6980  
info@polytec.co.jp

 **Polytec South-East Asia  
Pte Ltd**  
Blk 4010 Ang Mo Kio Ave 10  
#06-06 TechPlace 1  
Singapore 569626  
Tel. +65 64510886  
info@polytec-sea.com

 **Polytec China Ltd.**  
Room 1026, Hanwei Plaza  
No. 7 Guanghua Road  
Chaoyang District  
100004 Beijing  
Tel. +86 10 65682591  
info-cn@polytec.com