



A08615



Polarizer/Analyzer, A08615

A08712



Polarizer, A08712

*Dimensions in () are in mm

^Polarizer ring only. User must supply analyzer.

Ringlight Polarizers & Analyzers

The polarizing effect enhances contrast on highly reflective surfaces for Microscopy and Machine Vision Applications.

Standard Features


- The polarizer/analyzer reduces reflection and glare.
- The A08615, A08632, and A08652 polarizers were designed to easily snap-fit onto the ringlight. The analyzer is dropped through the top of the ringlight and positioned on top of the polarizer.
- To prevent the polarizer/analyzer from moving due to vibration, a knurled thumb-screw locks the polarizer assembly into position. In addition, it makes the holder easier to adjust. To further prevent vibration, a set-screw stops the polarizer disk from turning within the assembly. These features are found on polarizer / analyzer model numbers A08615 and A08627.

Polarizers & Analyzers


Part No.	Corresponding Ringlights
A08615	A08600, A08603, A08604, A08620, A08625, A08680, A22000, A22025
A0862	A08625 & Olympus Adapter, part number A08626
A08632	A08630, A08635
A08652*	A08650, A08660, A08670
A08705^	A08700
A08712^	A08710, A08715
A22042	A22040
A22052	A22050, A22060

Ringlight Polarizers & Analyzers


Dimensions in () are in mm




A08615
Polarizer/Analyzer




A08632
Polarizer/Analyzer




A08705
Polarizer/Analyzer




A22042
Polarizer/Analyzer




A08652
Polarizer/Analyzer



A08712
Polarizer



A22052
Polarizer/Analyzer



SCHOTT GLAS

Fiber Optics Division
P.O. Box 2480
55014 Mainz
Germany
Phone: +49(0) 61 31 / 66-0
Fax: +49(0) 61 31 / 66-7705
E-mail: fiber.optics@schott.com
www.schott.com/fiberoptics

SCHOTT FOSTEC

62 Columbus Street
Auburn, NY 13021
USA
Phone: +1 315 / 255-2791
Fax: +1 315 / 255-2695
E-mail: schott-fostec@us.schott.com
www.us.schott.com/fiberoptics
A Division of Schott North America, Inc.

